



PT.AVERON

— GLOBAL INDONESIA —

IMPORTER & SUPPLIER OF MEDICAL INSTRUMENTS



Stel Tec Care & Stroage

 www.averonindonesia.com




STEL TEC CARE & STORAGE








Storage Systems	Storage Rack Set	Item No.
	<p>For MIS instruments, standard For instruments with shaft length up to 37 cm PL963R storage rack perforated tray and</p>	ST 786
Parts	Description	Item No.
	Storage rack, standard dimensions: 487 x 250 x 153 mm	ST 787
	Tray, standard	ST 788
	Silicone mesh mat 240 mm x 250 mm	ST 789
	Silicone component	ST 790
	Trocar storage set universal silicone storage racks for trocars	ST 791
Container for Storage, Standard	Description	Item No.
	Container body external dimensions 592 x 274 x 187 mm	ST 792
	Lid for ST 777 external dimensions 588 x 285 x 36 mm	ST 793

Storage Systems	Storage Rack Set	Item No.
	<p>For MIS instruments, long For instruments with shaft length up to 42 cm, comprising: PL964R storage rack, long JF232R perforated tray and JF938 silicone mesh mat</p>	ST 794
Parts	Description	Item No.
	Storage rack, long dimensions: 587 x 250 x 153 mm	ST 795
	Tray, long	ST 796
	Silicone mesh mat 240 mm x 250 mm	ST 797
	Silicone component	ST 798
	Trocar storage set universal silicone storage racks for trocars	ST 799
Container for Storage, Extra Long	Description	Item No.
	Container body with perforation base external dimensions: 711 x 274 x 187 mm	ST 800
	Lid for ST 888 external dimensions: 708 x 280 x 35 mm	ST 801

Endoscopy Rack Storage System	Description	Item No.
	Endorack storage system, includes endorack, silicone mat, and 2"-high perforated basket.	ST 802
Container Accessories	Description	Item No.
	Endorack storage rack, holds 14 endoscopic instruments. Fits into 8" deep SterilContainer™ System.	ST 803
	Perforated basket, 21¼" long x 10" wide x 2¼" high	ST 804
	Silicone mat	ST 805
Container for Storage, Standard	Description	Item No.
	Container body external dimensions 592 x 274 x 187 mm	ST 806
	Lid for ST 999 external dimensions 588 x 285 x 36 mm	

Instrument Sterilization Trays & Containers	Description	Item No.
	Extra long mini container lid, for PreVac Steam and EtO	ST 807
	Extra long mini container lid, for PreVac Steam, EtO and low temp*	
	Extra long mini container bottom, for PreVac Steam and EtO	ST 808
	Extra long mini container bottom, for PreVac Steam, EtO and low temp*	
	Extra long mini basket with Lid, 20" long	ST 809

Endoscope Sterilization Trays & Containers	Description	Item No.
	Single scope tray, fits Full-Size SterilContainer™ System, 17" (45 cm) long	ST 810
	Double scope tray, fits Full-Size SterilContainer™ System, 17" (45 cm) long	ST 811
	Double scope tray with silicone supports, fits Full-Size SterilContainer™ System, 17" (45 cm) long	ST 812

Extra-Long Mini Size Container Accessories							
STS 25	STS 26	STS 27	STS 28	STS 29	STS 30	STS 31	STS 32
							
Silicone scope eyepiece holder	Silicone scope sleeve holder	Support for angled optics	Small silicone holder for 5 mm instruments	Small silicone holder for 7.5 mm instruments	Small silicone holder for 8.5 mm instruments	Small silicone holder for 10 mm instruments	Small silicone holder for 12 mm instruments



■ **Reliable...**

for patients and hospital staff. The Stellate SterilContainer System helps ensure proper storage, transportation and sterilization of surgical instruments.

■ **Cost-effective reusable container system...**

low operating expense and typical 10+ year service life.

■ **Easy to handle...**

simple, fast and efficient – lightweight aluminum and ergonomic handles.

■ **Clear identification match...**

using colored lids, indicator labels, identification labels and barcoding – clearly identifiable and standardized, easy to track and document.

In addition to a full line of containers, Stellate also offers storage systems, instrument racks and convenient disposable products to support all of your sterilization needs.

The Stellate Instrument Organization System is an advanced system of organizing, storing and protecting your instrument set during transportation and sterilization.

- Instruments organized with the Instrument Organization System are guarded against damage from instrument-on-instrument impact during transportation
- Instrument protection reduces the cost of repairs ^ set downtime
- Instrument protection increases instrument life ^ reduces replacement cost

In the O.R. a set organized with the Stellate Instrument Organization System can also help cases flow more efficiently by allowing staff to:

- Quickly identify and select the correct instruments
- Work directly from the organized set
- Replace used instruments into the correct slot making final instrument count easier

During sterilization, the Stellate Instrument Organization System allows:

- Better sterilant flow between instruments
- Consistent sterilization every time because your instruments are in the same place



CONTACT US

 +62 856-6381-381

 inquiry@averonindonesia.com

 sales@averonindonesia.com

 www.averonindonesia.com

 Axa tower - Kuningan City Lt 45-7A Jl Prof. DR. Satrio, Kuningan, Kecamatan Setiabudi, Kota Jakarta Selatan, Daerah Khusus Ibukota Jakarta 12940.