



Stel Tec Laparoscopic

 www.averonindonesia.com

STEL TEC LAPAROSCOPIC



Reusable 32 and 45 cm length laparoscopic instruments

- 32 and 45 cm lengths accommodate standard and advanced laparoscopic procedures, including Gynecology & Bariatrics.
- Rotating shafts with 360 degree rotation, one finger controls the tip placement.
- A wide variety of handle styles are available to accommodate a wide variety of surgeon preference.
- Extensive selection of tip styles.
- Offered in two styles to meet your cleaning requirements.
- Modular Laparoscopic instruments can be quickly disassembled for cleaning and visual inspection. Processing technicians spend less time removing bio-burden and the risk of HAI decreases when they can be certain that bio-burden is not trapped in the instrument. Easily reassembled for the surgeon prior to use in surgery

Single Piece with Flush-Port for quick flushing of bio-burden

Features:


- Durability and Longevity - Stainless steel construction, medical grade polymer coating guarantee long life, easy maintenance and safety.
- Precise & Ergonomic Handling - Jaw patterns are designed and machined for precision grasping, dissecting or cutting. Large-finger handles with tang and large diameter rotating drum with stainless steel roller bearings give you precise control and comfort while reducing hand fatigue.
- Optimal Cleaning & Maintenance - Flush ports and expanded push-rod clearance on non-modular instruments, minimum components and solid connections on modular instruments assures easy maintenance.
- Excellent Return on Investment - Completely reusable and fully-autoclavable.

STEL TEC LAPAROSCOPIC



	<u>Length</u>	<u>Style A:</u>	<u>Style B:</u>
5MM Ø, MIXTER GRASPER/ DISSECTOR RIGHT ANGLE (LONG)	33 CM	52255-07A	52255-07B
	36 CM	52236-07A	52236-07B
	45 CM	52155-07A	52155-07B
	24 CM	52224-07A	52224-07B
5MM Ø, MIXTER GRASPER/ DISSECTOR, SHORT RIGHT ANGLE	33 CM	52255-59A	52255-59B
	36 CM	52236-59A	52236-59B
	45 CM	52155-59A	52155-59B
	24 CM	52224-59A	52224-59B
5MM Ø, MIXTER GRASPER/ DISSECTOR 45 DEGREE ANGLE	33 CM	52255-08A	52255-08B
	36 CM	52236-08A	52236-08B
	45 CM	52155-08A	52155-08B
	24 CM	52224-08A	52224-08B
5MM Ø, DUCKBILL GRASPER/ DISSECTOR	33 CM	52255-09A	52255-09B
	36 CM	52236-09A	52236-09B
	45 CM	52155-09A	52155-09B
	24 CM	52224-09A	52224-09B
5MM Ø, NEEDLE-NOSE GRASPER/DISSECTOR	33 CM	52255-10A	52255-10B
	36 CM	52236-10A	52236-10B
	45 CM	52155-10A	52155-10B
	24 CM	52224-10A	52224-10B
5MM Ø, KELLY GRASPER/ DISSECTOR	33 CM	52255-11A	52255-11B
	36 CM	52236-11A	52236-11B
	45 CM	52155-11A	52155-11B
	24 CM	52224-11A	52224-11B
5MM Ø, LONG REACH, RIGHT ANGLE DISSECTOR, DBL- ACTION (REQUIRES LARGER CANNULA)	33 CM	52255-53A	52255-53B
	36 CM	52236-53A	52236-53B
	45 CM	52155-53A	52155-53B
	24 CM	52224-53A	52224-53B
5MM Ø, PETELIN DISSECTOR, 60 DEG. ANGLE LEFT, SINGLE-ACTION	33 CM	52255-12A	52255-12B
	36 CM	52236-12A	52236-12B
	45 CM	52155-12A	52155-12B
	24 CM	52224-12A	52224-12B
5MM Ø, PETELIN DISSECTOR, 30 DEG. ANGLE LEFT, SINGLE- ACTION	33 CM	52255-13A	52255-13B
	36 CM	52236-13A	52236-13B
	45 CM	52155-13A	52155-13B
	24 CM	52224-13A	52224-13B
5MM Ø, MICRO GRASPER/ DISSECTOR, FENESTRATED	33 CM	52255-14A	52255-14B
	36 CM	52236-14A	52236-14B
	45 CM	52155-14A	52155-14B
	24 CM	52224-14A	52224-14B

5mm Ø DISSECTORS & GRASPERS



	<u>Length</u>	<u>Style A:</u>	<u>Style B:</u>
5MM Ø, "WINER" STYLE FUNDUS MAXI GRASPER	33 CM	52255-326A	52255-326B
	36 CM	52236-326A	52236-326B
	45 CM	52155-326A	52155-326B
	24 CM	52224-326A	52224-326B
5MM Ø, "MCKERNAN" STYLE FUNDUS MAXI GRASPER	33 CM	52255-126A	52255-126B
	36 CM	52236-126A	52236-126B
	45 CM	52155-126A	52155-126B
	24 CM	52224-126A	52224-126B
5MM Ø, DOUBLE CUP FENESTRATED FUNDUS MAXI GRASPER	33 CM	52255-226A	52255-226B
	36 CM	52236-226A	52236-226B
	45 CM	52155-226A	52155-226B
	24 CM	52224-226A	52224-226B
5MM Ø, ENDO CLINCH FENESTRATED GRASPER	33 CM	52255-27A	52255-27B
	36 CM	52236-27A	52236-27B
	45 CM	52155-27A	52155-27B
	24 CM	52224-27A	52224-27B
5MM Ø, RETRO FACING RAT-TOOTH GRASPER w/ ATRAUMATIC CUP	33 CM	52255-28A	52255-28B
	36 CM	52236-28A	52236-28B
	45 CM	52155-28A	52155-28B
	24 CM	52224-28A	52224-28B
5MM Ø, ALLIGATOR GRASPER, 10MM JAW, RETRO FACING TEETH, DBL-ACTION	33 CM	52255-29A	52255-29B
	36 CM	52236-29A	52236-29B
	45 CM	52155-29A	52155-29B
	24 CM	52224-29A	52224-29B
5MM Ø, ALLIGATOR GRASPER, 20MM JAW, , RETRO FACING TEETH, DBL-ACTION	33 CM	52255-129A	52255-129B
	36 CM	52236-129A	52236-129B
	45 CM	52155-129A	52155-129B
	24 CM	52224-129A	52224-129B
5MM Ø, ALLIGATOR GRASPER, RETRO FACING TEETH, SINGLE-ACTION	33 CM	52255-30A	52255-30B
	36 CM	52236-30A	52236-30B
	45 CM	52155-30A	52155-30B
	24 CM	52224-30A	52224-30B
5MM Ø, CUP SHAPED GRASPER WITH RETRO FACING TEETH, DBL-ACTION	33 CM	52255-58A	52255-58B
	36 CM	52236-58A	52236-58B
	45 CM	52155-58A	52155-58B
	24 CM	52224-58A	52224-58B
5MM Ø, WAVE GRASPER, LONG BEVELED EDGE JAW, DBL-ACTION	33 CM	52255-32A	52255-32B
	36 CM	52236-32A	52236-32B
	45 CM	52155-32A	52155-32B
	24 CM	52224-32A	52224-32B

5mm Ø DISSECTORS & GRASPERS



	<u>Length</u>	<u>Style A:</u>	<u>Style B:</u>
5MM Ø, WAVE GRASPER, LONG BEVELED EDGE JAW, SINGLE ACTION	33 CM	52255-33A	52255-33B
	36 CM	52236-33A	52236-33B
	45 CM	52155-33A	52155-33B
	24 CM	52224-33A	52224-33B
5MM Ø, CLAW GRASPER, 2 X 3 TEETH, SINGLE-ACTION	33 CM	52255-36A	52255-36B
	36 CM	52236-36A	52236-36B
	45 CM	52155-36A	52155-36B
	24 CM	52224-36A	52224-36B
5MM Ø, CLAW GRASPER, 2 X 3 TEETH, W/ EXTRA ROW OF TEETH, DOUBLE-ACTION	33 CM	52255-136A	52255-136B
	36 CM	52236-136A	52236-136B
	45 CM	52155-136A	52155-136B
	24 CM	52224-136A	52224-136B
5MM Ø, RETRACTION GRASPER, LONG RIDGED JAWS, 2 X 3 ATRAUMATIC TEETH, DOUBLE-ACTION	33 CM	52255-34A	52255-34B
	36 CM	52236-34A	52236-34B
	45 CM	52155-34A	52155-34B
	24 CM	52224-34A	52224-34B
5MM Ø, RETRACTION GRASPER, LONG RIDGED JAWS, 2 X 3 ATRAUMATIC TEETH, SINGLE-ACTION	33 CM	52255-35A	52255-35B
	36 CM	52236-35A	52236-35B
	45 CM	52155-35A	52155-35B
	24 CM	52224-35A	52224-35B
5MM Ø, RETRACTION GRASPER, 2 X 3 ATRAUMATIC TEETH, DOUBLE-ACTION (JARIT STYLE)	33 CM	52255-134A	52255-134B
	36 CM	52236-134A	52236-134B
	45 CM	52155-134A	52155-134B
	24 CM	52224-134A	52224-134B
5MM Ø, RETRACTION GRASPER, 2 X 3 ATRAUMATIC TEETH, SINGLE-ACTION (JARIT STYLE)	33 CM	52255-135A	52255-135B
	36 CM	52236-135A	52236-135B
	45 CM	52155-135A	52155-135B
	24 CM	52224-135A	52224-135B
5MM Ø, ATRAUMATIC FENESTRATED GRASPER, CURVED, 16 MM DBL-ACTION JAW	33 CM	52255-37A	52255-37B
	36 CM	52236-37A	52236-37B
	45 CM	52155-37A	52155-37B
	24 CM	52224-37A	52224-37B
5MM Ø, FENESTRATED GRASPER, 16 MM DBL-ACTION JAW	33 CM	52255-38A	52255-38B
	36 CM	52236-38A	52236-38B
	45 CM	52155-38A	52155-38B
	24 CM	52224-38A	52224-38B
5MM Ø, FENESTRATED GRASPER, 20 MM DBL-ACTION JAW	33 CM	52255-39A	52255-39B
	36 CM	52236-39A	52236-39B
	45 CM	52155-39A	52155-39B
	24 CM	52224-39A	52224-39B

5mm BIOPSY FORCEPS & SCISSORS

	<u>Length</u>	<u>Style A:</u>	<u>Style B:</u>
	33 CM	52255-60A	52255-60B
	36 CM	52236-60A	52236-60B
	45 CM	52155-60A	52155-60B
	24 CM	52224-60A	52224-60B
5MM Ø, MICRO TIP SCISSORS, STRAIGHT, SINGLE-ACTION	33 CM	52255-61A	52255-61B
	36 CM	52236-61A	52236-61B
	45 CM	52155-61A	52155-61B
	24 CM	52224-61A	52224-61B
5MM Ø, MICRO TIP SCISSORS, CURVED, SINGLE-ACTION	33 CM	52255-62A	52255-62B
	36 CM	52236-62A	52236-62B
	45 CM	52155-62A	52155-62B
	24 CM	52224-62A	52224-62B
5MM Ø, METZENBAUM SCISSORS, STRAIGHT, LONG BLADE (20MM TIP)	33 CM	52255-63A	52255-63B
	36 CM	52236-63A	52236-63B
	45 CM	52155-63A	52155-63B
	24 CM	52224-63A	52224-62B
5MM Ø, METZENBAUM SCISSORS, STRAIGHT, SHORT BLADE (12MM TIP)	33 CM	52255-64A	52255-64B
	36 CM	52236-64A	52236-64B
	45 CM	52155-64A	52155-64B
	24 CM	52224-64A	52224-64B
5MM Ø, METZENBAUM SCISSORS, CURVED, LONG BLADE (18MM TIP)	33 CM	52255-65A	52255-65B
	36 CM	52236-65A	52236-65B
	45 CM	52155-65A	52155-65B
	24 CM	52224-65A	52224-65B
5MM Ø, METZENBAUM SCISSORS, HARDENED STEEL, CURVED BLADE (16MM TIP)	33 CM	52255-66A	52255-66B
	36 CM	52236-66A	52236-66B
	45 CM	52155-66A	52155-66B
	24 CM	52224-66A	52224-66B
5MM Ø, METZENBAUM SCISSORS, CURVED, SHORT BLADE (11MM TIP)	33 CM	52255-166A	52255-166B
	36 CM	52236-166A	52236-166B
	45 CM	52155-166A	52155-166B
	24 CM	52224-166A	52224-166B
5MM Ø, PERITONEAL SCISSORS, SINGLE-ACTION	33 CM	52255-67A	52255-67B
	36 CM	52236-67A	52236-67B
	45 CM	52155-67A	52155-67B
	24 CM	52224-67A	52224-67B
5MM Ø, HOOK SCISSORS, TS, EXTRA HARD STEEL	33 CM	52255-76A	52255-76B
	36 CM	52236-76A	52236-76B
	45 CM	52155-76A	52155-76B
	24 CM	52224-76A	52224-76B

5mm Ø DISSECTORS & GRASPERS

	<u>Length</u>	<u>Style A:</u>	<u>Style B:</u>
	33 CM	52102-02A	52102-02B
	36 CM	52103-02A	52103-02B
	45 CM	52101-02A	52101-02B
	33 CM	52102-06A	52102-06B
	36 CM	52103-06A	52103-06B
	45 CM	52101-06A	52101-06B
	33 CM	52102-07A	52102-07B
	36 CM	52103-07A	52103-07B
	45 CM	52101-07A	52101-07B
	33 CM	52102-59A	52102-59B
	36 CM	52103-59A	52103-59B
	45 CM	52101-59A	52101-59B
	33 CM	52102-08A	52102-08B
	36 CM	52103-08A	52103-08B
	45 CM	52101-08A	52101-08B
	33 CM	52102-16A	52102-16B
	36 CM	52103-16A	52103-16B
	45 CM	52101-16A	52101-16B
	33 CM	52102-17A	52102-17B
	36 CM	52103-17A	52103-17B
	45 CM	52101-17A	52101-17B
	33 CM	52102-21A	52102-21B
	36 CM	52103-21A	52103-21B
	45 CM	52101-21A	52101-21B
	33 CM	52102-22A	52102-22B
	36 CM	52103-22A	52103-22B
	45 CM	52101-22A	52101-22B
	33 CM	52102-23A	52102-23B
	36 CM	52103-23A	52103-23B
	45 CM	52101-23A	52101-23B

10mm FORCEPS, BIOPSY, TENACULUMS

	<u>Length</u>	<u>Style A:</u>	<u>Style B:</u>
	33 CM	52102-24A	52102-24B
	36 CM	52103-24A	52103-24B
	45 CM	52101-24A	52101-24B
	33 CM	52102-25A	52102-25B
	36 CM	52103-25A	52103-25B
	45 CM	52101-25A	52101-25B
	33 CM	52102-26A	52102-26B
	36 CM	52103-26A	52103-26B
	45 CM	52101-26A	52101-26B
	33 CM	52102-36A	52102-36B
	36 CM	52103-36A	52103-36B
	45 CM	52101-36A	52101-36B
	33 CM	52102-37A	52102-37B
	36 CM	52103-37A	52103-37B
	45 CM	52101-37A	52101-37B
	33 CM	52102-38A	52102-38B
	36 CM	52103-38A	52103-38B
	45 CM	52101-38A	52101-38B
	33 CM	52102-45A	52102-45B
	36 CM	52103-45A	52103-45B
	45 CM	52101-45A	52101-45B
	33 CM	52102-94A	52102-94B
	36 CM	52103-94A	52103-94B
	45 CM	52101-94A	52101-94B
	33 CM	52102-46A	52102-46B
	36 CM	52103-46A	52103-46B
	45 CM	52101-46A	52101-46B
	33 CM	52102-47A	52102-47B
	36 CM	52136-47A	52136-47B
	45 CM	52101-47A	52101-47B
	33 CM	52102-64A	52102-64B
	36 CM	52136-64A	52136-64B
	45 CM	52101-64A	52101-64B
	33 CM	52102-67A	52102-67B
	36 CM	52136-67A	52136-67B
	45 CM	52101-67A	52101-67B

STEL TEC LAPAROSCOPIC



TYPE 1:
INSULATED WITH RF-
POST



TYPE 2:
INSULATED W/ RF-
POST & RATCHET



TYPE 3:
NON-INSULATED
STAINLESS STEEL



TYPE 4:
NON-INSULATED
STAINLESS STEEL W/
RATCHET



TYPE 5:
LIGHTWEIGHT
CARBON-FIBER W/ RF
POST



TYPE 6:
LIGHTWEIGHT
CARBON-FIBER
w/ RF-POST &
RATCHET



TYPE 7:
LIGHTWEIGHT
CARBON-FIBER
W/ RF POST,
RATCHET /
RATCHET
DEFEATING



TYPE 8:
LIGHTWEIGHT CARBON
-FIBER W/ RF POST,
FRENCH ON/OFF
RATCHET



TYPE 9:
LIGHTWEIGHT
CARBON-FIBER W/ RF
POST, W/ LARGE
FINGER LOOPS



TYPE 10:
LIGHTWEIGHT
CARBON-FIBER W/ RF
POST, RATCHET /
RATCHET DEFEATING
W/ LARGE FINGER
LOOPS



TYPE 11:
STAINLESS STEEL RING
AXIAL HANDLE



TYPE 12:
STAINLESS STEEL RING AXIAL
HANDLE w/ RATCHET



TYPE 13:
STAINLESS STEEL AXIAL HANDLE
w/ PRECISE ENGAGEMENT
RATCHET



TYPE 14:
INSULATED STAINLESS
STEEL LINVATEC-WECK
STYLE W/ LONG
SHANKS



TYPE 15:
NON-INSULATED
STAINLESS STEEL
LINVATEC-WECK STYLE
W/ QUICK RELEASE
RATCHET & LONG
SHANKS



TYPE 16:
NON-INSULATED
STAINLESS STEEL W/
SPRING ASSIST
CLOSING



TYPE 18:
CARBON-FIBER W/ RF
POST



TYPE 18:
CARBON-FIBER W/ RF
POST, RATCHET /
RATCHET DEFEATING



CONTACT US

 +62 856-6381-381

 inquiry@averonindonesia.com

 sales@averonindonesia.com

 www.averonindonesia.com

 Axa tower - Kuningan City Lt 45-7A Jl Prof. DR. Satrio, Kuningan, Kecamatan Setiabudi, Kota Jakarta Selatan, Daerah Khusus Ibukota Jakarta 12940.